



Republic of the Philippines
Department of Agriculture
BUREAU OF FISHERIES AND AQUATIC RESOURCES
REGIONAL FISHERIES LABORATORY XII
J. Catolico St., Lagao, General Santos City

Document Type Standard Operating Procedures	Revision No: 0	Date Adopted: 13 Aug 2019	Page No: Page 1 of 2
Document Code: SOP-02-04	WASHING GLASSWARE		

Introduction

This document provides instructions for washing glassware which have been used by or are for use in analysis in the BFAR RFL XII. It ensures a uniform cleaning process to minimize, if not completely eliminate, contaminants.

1. Responsibilities

All Laboratory Staff

2. Procedure

- 2.1 Rinse glassware thoroughly with running water.
- 2.2 Soak glassware in laboratory detergent soaking solution for at least 10-30 minutes for those used in chemistry laboratory and overnight for those used in microbiology and molecular diagnosis laboratory.
- 2.3 Rinse detergent off thoroughly with running water.
- 2.4 Rinse thoroughly and repeatedly with tap water.
- 2.5 Acid wash if necessary.
- 2.6 Rinse thoroughly with distilled water.
- 2.7 Let glassware dry before storing them in their respective shelves or cupboards immediately after. Storage areas are marked visibly for easy access.

3. Materials

- 3.1 Laboratory detergent
- 3.2 Laboratory detergent soaking solution
- 3.3 Disinfectant cleaning solution
- 3.4 Acid Solution
- 3.5 Basins and pails for soaking
- 3.6 Distilled water
- 3.7 Wash bottles

Prepared by: MA. MAURICE A. BANQUILING Quality Assurance Manager	Reviewed and Approved by: EUGENE GAY B. JAMORA Laboratory Manager
--	---



Republic of the Philippines
Department of Agriculture
BUREAU OF FISHERIES AND AQUATIC RESOURCES
REGIONAL FISHERIES LABORATORY XII
J. Catolico St., Lagao, General Santos City

Document Type Standard Operating Procedures	Revision No: 0	Date Adopted: 13 Aug 2019	Page No: Page 2 of 2
Document Code: SOP-02-04	WASHING GLASSWARE		

4. Frequency

- 4.1 These instructions must be followed every time glassware are cleaned.
- 4.2 The Laboratory Analysts are responsible for the implementation of this procedure.
- 4.3 The Technical Manager ensures that the objectives of this procedure are met.

5. Document History

Revision No.	Issue Date	Original Entry	Revised as	Recall Date

Prepared by:

MA. MAURICE A. BANQUILING
Quality Assurance Manager

Reviewed and Approved by:

EUGENE GAY B. JAMORA
Laboratory Manager