



Republic of the Philippines
Department of Agriculture
BUREAU OF FISHERIES AND AQUATIC RESOURCES
REGIONAL FISHERIES LABORATORY XII
J. Catolico St., Lagao, General Santos City

Document Type Work Instruction	Revision No: 0	Date Adopted: 13 Aug 2019	Page No: Page 1 of 2
Document Code: WI-01-MIC-01	SAMPLE PREPARATION		

Introduction

This document provides a standard procedure to ensure laboratory samples taken from the whole sample submitted are of homogeneous quality. This procedure is applicable to fish and fishery products submitted for Microbiological Analysis to the Regional Fisheries Laboratory.

1. Responsibilities

- Laboratory Analyst
- Technical Manager
- Quality Assurance Manager
- Customer Service Officer

2. Procedure

- As soon as the sample arrives at the laboratory, the analyst should note its general physical condition (sealed).
- A frozen sample can be thawed at 2-5°C within 18 h. If rapid thawing is desired, thaw the sample at less than 45°C for not more than 15 min.
- Weigh 25 g sample into a sterile stomacher bag.
- Proceed as indicated in the specific tests to be performed.

Prepared by: CHARISSE ZIANNE N. LIBRES Laboratory Analyst	Reviewed by: MA. MAURICE A. BANQUILING Quality Assurance Manager	Approved by: EUGENE GAY B. JAMORA Laboratory Manager
---	--	--



Republic of the Philippines
Department of Agriculture
BUREAU OF FISHERIES AND AQUATIC RESOURCES
REGIONAL FISHERIES LABORATORY XII
J. Catolico St., Lagao, General Santos City

Document Type Work Instruction	Revision No: 0	Date Adopted: 13 Aug 2019	Page No: Page 2 of 2
Document Code: WI-01-MIC-01	SAMPLE PREPARATION		

3. Document History

Revision No.	Issue Date	Original Entry	Revised as	Recall Date

Prepared by: CHARISSE ZIANNE N. LIBRES Laboratory Analyst	Reviewed by: MA. MAURICE A. BANQUILING Quality Assurance Manager	Approved by: EUGENE GAY B. JAMORA Laboratory Manager
---	--	--